



New runners of ophthalmic market
BIO OPTICS

New runners of ophthalmic market

BIO OPTICS

Product Catalog



No.1215, SK V1 Tower, Galmachi-Str 288-14, 13201,
Jungwon-gu, Seongnam-si, Gyeonggi-do, Republic of Korea
Tel. +82 31 747 4850 / Fax. +82 31 747 4820
E-mail : Zevany1003@gmail.com / Totour1187@naver.com
www.biooptics.co/en



RETINALAB
Medical Company



Distributor in Azerbaijan

Hello, eyes **BIO Plug™**

Dissolvable Lacrimal Punctum Plug



Hello, eyes **BIO Intubation Set™**

Lacrimal Intubation Set



5 Sets / Box

Features

- New state art of lacrimal intubation set
- Bio compatible silicon tubing
- Patent acquired for South Korea
- Made of one of best silicon material from German company

REF Num. BIO Plug 0.3 / BIO Plug 0.4 / BIO Plug 0.5




10 Pcs / Box

Features


- Made of Polydioxanone as material
- Violet Color
- Never Touch, and less irritate the eye
- Comfortable after proper placement
- Absorption of plug in human body approximately 180 days
- Available in 0.3mm, 0.4mm and 0.5mm diameters.

Indications

- Dry eye syndrome
- Refractive surgery
- Contact lens wear
- Seasonal allergy
- Cataract surgery



Stainless steel probes with silicone tubing




BS-1001 [23G] / BS-1002[27G]

Stainless Probe with Olive Tip

11cm length x 0.65 mm diameter [BS-1001] / 11cm length x 0.45 mm diameter [BS-1002]

Silicon Tubing

30cm length x 0.94 mm diameter [BS-1001] / 30cm length x 0.64 mm diameter [BS-1002]




BS-2001 [23G] / BS-2002 [27G]

Stainless Probe without Olive Tip

11cm length x 0.65 mm diameter [BS-2001] / 11cm length x 0.45 mm diameter [BS-2002]

Silicon Tubing

30cm length x 0.94 mm diameter [BS-2001] / 30cm length x 0.64 mm diameter [BS-2002]



BS-3001

Malleable stainless steel and nylon sheath.

The rounded tip, nylon sheath diameter offers a smooth transition to silicon tubing

Stainless steel probe covered in silicon tube - 15 cm x 0.51 mm diameter


Silicon Tube with knotted with nylon sheath - 33 cm x 0.64 mm diameter.



BS-4001 Retrieval device

Stainless steel intubation retrieval device with hook shape head


Hello, eyes **BIO Ophthalmic Surgical Knife**



New feature of Laser marking depth indicator in keratome knives.


Knife products introduction.

1. BIO Keratome Knife




Model No.	Length (mm)	Blade Length (mm)	Blade Diameter (mm)	Description
B36-2200	1350	200	2.20	Phaco-Slit Knife, Angled, Double-Bevel
B36-2250	1350	200	2.25	
B36-2400	1350	200	2.40	
B36-2600	1350	200	2.60	
B36-2750	1350	200	2.75	
B36-2800	1350	200	2.80	
B36-3000	1350	200	3.00	

2. BIO MVR Knife



Model No.	Length (mm)	Blade Length (mm)	Blade Diameter (mm)	Description
B36-1900	1350	300	0.75	19G MVR Knife, Straight
B36-1901	1350	300	0.75	19G MVR Knife, Angled
B36-2000	1350	300	0.55	20G MVR Knife, Straight
B36-2001	1350	300	0.55	20G MVR Knife, Angled

3. BIO Scleral Knife




Model No.	Length (mm)	Blade Length (mm)	Blade Diameter (mm)	Description
B36-5700	1350	200	3.2	200 / 3.2

4. BIO Side Port Knife



Model No.	Length (mm)	Blade Length (mm)	Blade Diameter (mm)	Description
B36-7511	1350	200	1.0	15° stab knife, Straight
B36-7513	1350	200	1.5	
B36-7515	1350	200	3.0	

5. BIO Crescent Knife



Model No.	Length (mm)	Blade Length (mm)	Blade Diameter (mm)	Description
B36-4805	1350	200	2.0	Bevel-up Crescent Knife, Straight
B36-4806	1350	200	2.0	Bevel-up Crescent Knife, Angled

Products

Hello, eyes BIO Color Tear Test



Features

- Easy to use
- Blue-black dye line
- Different shape for right and left eyes
- Made of the best whatman filter paper
- paper from German

Critical values

- Normal which is 10mm-30mm
- wetting of the paper after 5 Min
- Dry eye patient which is <10 mm
- wetting of the paper after 5 Min

Hello, eyes BIO Micro Syringe Filter



Indications

- Refractive Surgery
- Cataract Surgery
- Retina, Glaucoma and other Ophthalmic surgery

Specifications

- 0.2mm x 25mm
- Safety to PH 1~14
- Hydrate membrane filter
- Gamma Radiation

Hello, eyes EYE CLEAN LID CLEANER



Features

- The first tea tree oil extract product in South Korea
- Excellent effect to relieve the dry eye syndrome
- Unique formulation for eliminating Demodex mites by tea tree oil extract
- Contain no Alcohol, Fragrance, Preservative materials and no paraben.
- There is no needs to rinse or remove after use for the comfort for the customer
- Eyeclean eyelids cleaner made of 100 % cotton patch
- Never irritate delicate eyelids and eyelashes

Hello, eyes EYE CLEAN CLEAN DEX



Features

- Different type of Eyeclean eyelids cleaner
- Highly recommended for the customer, who wants to use ling period of eye cleaner
- It is like form type so pull the trigger and use the inner liquid with enclosed cotton pads
- Contains 100 % cotton pads.

Indications

- After Make-up removal, complete wipe out of oil residue of eye Make-up.
- Relieve the dry eye syndrome
- Inflamed and itchy eye
- Helps to eliminate the Demodex mites
- Relieve Blepharitis, Stye such ophthalmic diseases.
- Comforting for contact lens
- Before and After cataract and other ophthalmic surgeries.
- Cleaning eyelids and eyelashes.