


BASE

series

PROCEDURE CHAIR

This procedure chair can be used as ophthalmic table, eye examination table and so on. It's a multifunctional physiotherapy chair and suitable for a wide range of treatment & other therapy.



Scan the QR code
For product application details

4 Motors

Warranty

Wear Resistance

Certificate

250kg/551lbs

More



Product features

- Overall Dimensions: Length 198cm/78", Width 60cm/23.6"
- Height Adjustment: From 59cm/23.2" to 89cm/38.6"
- Electric adjustable height, back section, leg section, seat section by handset operation.
 - Pneumatic adjustable head section
- Head section: -15° to +35° by gas spring controlled
- Back section: -5° to +85° by motor controlled, leg section: 0° to +85° by motor controlled
- Seat section: 0° to +35° by motor controlled
- Contoured comfortable table top and flip up armrests, removable padded wing armrests
- Debris tray included
- Unobstructed base design. All-welded base and top frame
- Motors: Reliable column motor.
- Upholstery Foam: 50mm/2" thickness with a density of 36kg/m³ (2.25lbs/ft³)
- Vinyl: Fire-retardant PVC vinyl (BS 5852 compliant, 150,000 abrasion cycles)
- Castor System: One lever controls four retractable castors (Ø125mm/5")
- Mains Supply: 100-240V AC, 50/60 Hz, Class II
- Weight Capacity: 250kg/551lbs, optional: 280kg/617lbs
- Warranty: 5 years
- Certifications: ISO13485, CE, UKCA, TGA, FDA approved

Technical data

Padding thickness
·50mm / 2"

Foam density
·36kg/m³ / 2.25lbs/ft³

Castors
·φ125mm / φ5"

Product weight
·111kg / 245lbs

Packing size
·183x82x49cm / 72.1"x32.3"x19.3"

Optional accessories

Padded side rails

Paper roll holder

Disposable paper

IV pole

Safety lock

Battery

Product model

BASE SERIES

Product Name

ELT35PRO

Control mode	Motor
Height adjustment	Handset
Lifting/Down speed	30s
LxWxH	198x60x(59-89)cm 78"x23.6"x(23.2"-38.6")

Colour of vinyl

Sky blue

Navy blue

Charcoal grey

White

Cement grey

Black

Colour of frame

RAL 9003

RAL 9002

RAL 7024